# CFC - Laborants en chimie

## Anglais 1999

### Contexte

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## 47. Effect of Temperature Change on Equilibrium: Cobalt Complex

An equilibrium system involving the dehydrated-hydrated cobalt complex is produced. When this system is heated, a color change from pink to blue indicates a shift of equilibrium to the right. When the solution is cooled, the color change from blue to pink indicates a shift to the left.


### Procedure

Wear safety goggles and disposable plastic gloves.  $\text{CoCl}_2$  is toxic.

1. Place 100 mL of  $\text{CoCl}_2$  solution in a 250-mL beaker.
2. Add concentrated HCl until the solution changes from pink to blue. Use a face shield and gloves when you use concentrated HCl.
3. Divide the solution into three smaller beakers and treat them as follows:
  - a. Place one beaker on a hot plate.
  - b. Place one beaker in an ice bath.
  - c. Leave one beaker at room temperature as a standard.
4. After a few minutes, show that the heated sample has turned a darker blue and that the cooled sample has turned a light pink.

### Reactions

1. This reaction involves the following equilibrium:

2. Addition of heat causes a shift of equilibrium toward products, the blue solution.
3. Cooling causes a shift of equilibrium to the left, the pink hydrated complex.

### Solutions

1. The  $\text{CoCl}_2$  solution is 0.4 M: Dissolve 5.2 g per 100 mL of water.
2. The HCl solution is concentrated.

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a) Why do you wear plastic gloves, a face shield and safety goggles ?

b) Is it dangerous to work in a lab ?

c) Is it difficult to realize such an experience ?

d) Do you always work with a white coat ?

e) Would you recommend me to learn this job ?

f) How long does it take to learn this job ?

g) Did you enjoy your apprenticeship ?